


How to Collect Your Nasal Swab Sample

Thank you for getting a COVID-19 diagnostic test today! Please follow the steps below to collect your nasal swab sample. We are here to help if you have any questions or need help.


1 Open Nasal Swab

Remove the nasal swab from the wrapper by pulling the two ends apart (like you would open a band-aid). Be careful to only touch the handle, not the tip.




2 Swab Your Nose

Gently put the soft tip of the swab into one nostril, until you feel you can't push it further. Rub the swab around 4 times. Take out the swab slowly. Repeat this step in your other nostril with the same swab.



3 Put Swab in the Tube

Place the swab in the tube we gave you, tip first. When the tip touches the bottom, bend the swab handle to break it off. Give your tube to a member of our staff, who will share next steps.



We will send your sample for lab testing. When your results are ready, you will receive them via MyChart.